

Rapid Strain and Early Process Development for Recombinant Erwinia Asparaginase (Rylaze™)

Ligand Pharmaceuticals, San Diego, California

Torben Bruck, Russell Coleman, Jason Payne, Keith Haney, Diane Retallack

Rapid and comprehensive strain selection and an early process development program utilizing the Pelican Expression Technology™ were key to establishing the foundation for late-stage success for Rylaze™, a recombinant *Erwinia chrysanthemi* asparaginase for the treatment of acute lymphoblastic leukemia (ALL) or lymphoblastic lymphoma (LBL) in patients who have developed hypersensitivity to *E. coli*-derived asparaginase. Pelican Expression Technology™, a *Pseudomonas fluorescens*-based protein expression system is a robust and scalable platform for recombinant protein production.


Rapid Strain Construction, Expression, and Analysis



- 2 Codon Optimized Genes
- 2 Plasmids
- 4 Transcription Promoters
- 3 RBS/ Translation Initiation Strengths
- 53 P. Fluorescens Native Signal Peptides
- 60 Folding Helper Overexpression Host Strains
- 150 Protease Deletion Host Strains

- Toolbox components combined to produce **thousands of unique expression strains** constructed and screened in parallel using automated workflows
- Interactions between toolbox components are explored to quickly identify production strain candidates using high throughput analysis to assess target quantity and quality
- Parallel automated processing enables rapid strain screening in 9 weeks

ALL/LBL & Role of Asparaginases in Treatment



Acute Lymphoblastic Leukemia (ALL) is the most common form of childhood cancer affecting nearly 6000 patients per year in the US but also among the most curable!

Lymphoblastic Lymphoma (LBL) is a rare, fast-growing, aggressive sub-type of Non-Hodgkin's lymphoma most often seen in teenagers and young adults


Asparaginases are a key component of multi-agent chemotherapeutic regimen for the treatment of ALL and LBL and full and effective asparaginase administration is essential to best outcomes in patients

Until recently and due to world-wide supply challenges, patients with hypersensitivity to an *E. coli*-derived asparaginase had a critical unmet need for an immunologically distinct asparaginase product that would allow them to continue and complete their optimal treatment plan

Tetrameric asparaginase *Erwinia chrysanthemi*

1. National Cancer Institute, <https://www.cancer.gov/types/childhood-cancers/child-adolescent-cancers-fact-sheet>



Unique high titer crisantaspase expression strains identified at 96-well scale




- ▶ HTP 96-well strain screening identified strains producing active crisantaspase (or L-asparaginase) greater than 1 g/L soluble monomer was observed
- ▶ Construction of *P. fluorescens* host strains deficient in native L-asparaginase completed (chromosomal KOs)
- ▶ 5 total strains, representing a diversity of expression strategies, scaled to 2L fermentations for evaluation under multiple induction conditions
 - 4 periplasmic and 1 cytoplasmic localization strains

Fermentation Assessment

- ▶ Cytoplasmically localized PF743 (Strain #4) showed highest overall soluble titers of monomer ranging from 19 to 23 mg/L
- ▶ PF743 Strain #4 was cultivated under 8 different induction conditions (factors: wet cell weight at induction, IPTG concentration, pH, temperature)




High Quality Material for Preclinical Development




- Methods used for evaluation:
- I. Activity using Nessler-based method
 - II. Tetramer by SEC
 - III. Western/SDS-PAGE
 - IV. Expected mass by LC-MS
 - V. Purity by RP-HPLC
 - VI. CD/IF

LC-MS: Intact Mass



- ✓ Small-scale expression and purification of crisantaspase using the Pfenex Expression Technology successfully demonstrated
- ✓ Greater than 2.5 grams of *P. fluorescens* produced crisantaspase meeting initial purity and potency targets was produced
- ✓ Chosen production strain moved forward into process development


Upstream Process Development



- ▶ 5 strains advanced
- ▶ Up to 9 variable induction conditions screened using DoE approach
 - Typical factors tested
 - Temp.
 - pH
 - Inducer level
 - Wet cell weight
- ▶ Selection of robust production strain for further development

134 x 2L fermentations in 6 Design-of-experiment rounds to optimize conditions

YIELD RESPONSE: Surface plots of key variables



Phase 3 ready process enabled transfer directly to the intended commercial scale

28 x 20L fermentations to optimize conditions and support Downstream Processing Development